

About Esco
World-Class Test Facility 02
American Standard ASHRAE 110-1995
European Standard EN14175 03
Energy Efficiency
Benefits For All Stakeholders 05
Frontier Acela _® High Performance
Fume Hood Laminar Topography and
Perfect Protection Perimeter 06
Computational Fluid Dynamics 07
Safety Certified and Tested 08
General Specifications
Guide to Models
Technical Specifications 10
Options and Accessories 11
Frontier Acela Base Cabinetry (EBA)12
Esco Resinate

Ordering a Frontier Acela Fume Hood ... 15





About Esco

Since Esco was founded in 1978 our company has earned a reputation for innovation in the worldwide laboratory equipment and pharmaceutical equipment industry. Today, Esco has emerged as a market leader in containment, clean air, pharmaceutical, and laboratory equipment technologies with active sales in more than 100 countries and direct company offices in the top ten geospecific markets.

From our headquarters in Singapore, Esco directs a highly efficient research, product development, manufacturing and customer service program. We are the only company in our market that is completely configured to export most of what we manufacture.

Our many languages and cultures, customs and traditions, and modern business management techniques blend into a single effort focusing on customer service, one customer at a time. As you learn more about Esco, you will understand why World Class. Worldwide. is more than a phrase. It's part of who we are, where we are from and where we are going.

World-Class Test Facility

Esco is proud to be one of the few manufaturers in the world with a test facility capable of testing hoods to both ASHRAE 110-1995 (US) and the EN14175-3 (European) standards. Esco's Fume Hood Test Laboratory was designed with the assistance of Tintschl Engineering AG, a specialist consulting firm from Germany. It has controlled relative humidity, room temperature and pressure for optimum test conditions. Esco is also one of the few companies to routinely sample and subject production fume hoods to a battery of containment and safety tests. All custom fume hoods with modified dimensions are also tested in our laboratory to ensure containment before delivery.

A Complete Range

Frontier Acela_® is the world's first high performance / low flow fume hood product range with an ergonomically angled front. Frontier Acela is available in Benchtop, Perchloric Acid, Radioisotope, Acid Digestion and Floor Mounted configurations, in a wide range of widths, with a full range of options and accessories.

The ASHRAE 110-1995 standard is a comprehensive method for evaluating the operator safety of fume hoods by determining quantitatively and repeatably how well the fume hood contains vapors released in the workzone. First published in 1985 and extensively revised in 1995, this standard employs a set of rigorous tests to evaluate hood performance:


1. Airflow Visualization

 Local and gross airflow visualization tests conducted to observe airflow patterns inside work chamber

Grading	Observations
Fail	Smoke observed escaping from the hood
Poor	Reverse flow of smoke near opening Lazy flow into opening along boundary Observed potential for escape
Fair	Some reverse flow, not necessarily at opening No visible escape
Good	No reverse flows Active flow streams into hood around boundary

2. Face Velocity Measurements

Face velocity and uniformity tested at 100% and design sash opening position



3. Tracer Gas Containment

- Tracer gas containment test: SF6 released at 4L/min inside work chamber
- Consists of 3 parts
 - Static gas leak conducted at 3 positions (left, right, centre)
 - Face hood surface scan
 - Sash movement effect
- Test criteria: leakage of less than 0.05ppm as-manufactured




European Standard EN 14175

EN 14175 is a harmonized European standard which supercedes the former national standards of Germany, the UK and France. A key element of the EN 14175, which is not present in the American ASHRAE 110 standard, is the robustness test, which simulates airflow disturbance in front of the hood.

1. Inner Grid Test

 A set of 9 sampling tubes are arranged in a 3 x 3 grid, at a distance of 30 mm from the measurement



plane

• Leakage will be measured at 6 different grid positions

2. Outer Grid Test


- Outer grid testing measures the leakage of tracer gas across the entire sash opening
- Sash is closed then opened during the test to investigate dynamic effect of sash movement



3. Robustness Test

- For the robustness test, a 0.4 x 1.9 m plate crosses the front of the fume hood at 1 m/s for 6 times
- Provides an indication of the sensitivity of the fume hood towards external disturbances







3

Fume hoods — essential safety devices used in laboratory environments — are highly energy-intensive, each one consuming more energy than three homes in an average U.S. climate. Depending on climate and system design, estimated energy costs for fume hoods range up to **US\$9000 annually¹**, based on face velocities of 0.5m/s (100fpm) at full sash open position for a 1.8m (72.0") hood.

Variable Air Volume (VAV) is one of the various approaches presently employed to reduce hood energy consumption. The table below compares conventional hoods, VAV hoods, and the Esco Frontier Acela_® High Performance Low Flow Hood.

	Conventional Fume Hood	Variable Air Volume (VAV) Fume Hood	High Performance Low Flow Fume Hood	
Working Principle	0.5 m/s (100 fpm) @ full open sash position	0.5 m/s (100 fpm) @ all sash positions with sophisticated control system	0.3 m/s (60 fpm) @ 457 mm (18") sash opening using advanced aerodynamic designs	
Initial Cost	Low	High	Medium	
Running Cost	Very High	Medium (VAV Maintenance)	Low	
Ease of Installation, Commissioning and Maintenance	Easy	Difficult	Easy	


¹ Energy use and savings potential for laboratory fume hoods, Evan Mills, Dale Sartor; Energy, 2003



Compared with conventional hoods, Esco Frontier Acela_® operates safely at 0.3 m/s (60 fpm) at 457 mm (18.0") or full open sash position while maintaining excellent ASHRAE and EN containment. Exhaust volume reductions of up to 58% may be


achieved without compromising safety. **This translates into an annual operating cost savings of up to US\$5600**. Unlike VAV systems the Esco Frontier Acela_® is easy and inexpensive to install, commission and maintain.

	Exhaust			
Fume Hood Width	Frontier Acela _® 0.3 m/s (60 fpm) @ 457 mm (18")	Conventional Fume Hood 0.5 m/s (100 fpm) @ full sash open	% Reduction in Exhaust Volume	
1.2m (4')	541 cmh (316 cfm)	1296 cmh (763 cfm)	58%	
1.5m (5')	777 cmh (457 cfm)	1400 cmh (824 cfm)	44%	
1.8m (6')	872 cmh (510 cfm)	1543 cmh (908 cfm)	43%	
2.4m (8')	1203 cmh (708 cfm)	2372 cmh (1396 cfm)	49%	





Benefits For All Stakeholders



Health and Safety Authorities

- ASHRAE 110 and EN 14175 tested at face velocities down to 0.3m/s (60fpm)
- UL1805 Listed
- Sash-lock and creep-down mechanisms ensure operator safety
- Optional Sentinel XL Airflow Monitor

Facility Owners




- Operates at 0.3 m/s (60 fpm) at 457 mm (18.0") operating height
- Energy savings of up to US\$5600 annually for each fume hood



Lab Personnel

- Certified safety
- Precisely-tuned aerodynamics
- High sight line
- Angled front
- Excellent ergonomics



Lab Designers and Architects


- Neutral color scheme
- Futuristic design
- Detailed architectural specifications and drawings



Installation and Service Personnel

- Robust tri-wall side pan construction
- Factory-fitted service fixtures are pre-plumbed
- Chain and sprocket sash
- Wide range of field-installed accessories








Ergonomic Sash Handle


Gently directs air into the hood without sacrificing visibility.






Tapered Fiberglass Exhaust Collar

Transitions exhaust air quietly and enhance airflow uniformity.




Functionally Robust Bypass

The unique design provides a robust stream of bypass air into the hood cavity.



Hot Zone Baffles

The unique Hot Zone Baffles draw most contaminants back in single pass displacement of the air. Thermal Heat Relief is quickly achieved.



Aerodynamic Foil Entry


Provides maximum airflow "sweep" on the critical boundary layer near to the work surface level to reduce turbulence and eliminate backflow.

Computational Fluid Dynamics


Computational Fluid Dynamics (CFD) modelling is employed in the development of Esco clean air and containment devices. Laminar Topography™ on Frontier Acela® Fume Hoods was developed with computational fluid dynamics modelling in the Esco Research and Development Center. The main thrusts of the project were improved airflow uniformity, enhanced safety, reduction in noise levels, and energy consumption.

First, engineering teams conceptualized possible designs, and, instead of building physical models, utilized CFD to simulate airflow patterns, pressurizations and visualize possible areas of turbulence. This allowed a large number of iterations of the airfoil, sash handle, baffle, bypass and exhaust collar to be evaluated. Finally, physical prototypes were constructed, tested, and the best design combination selected for production.


CFD has allowed us to effectively reduce the vortex in conventional fume hood designs to the minimum, resulting in a safe yet energy-saving fume hood design.



Conventional Fume Hood



Frontier Acela_® High Performance Low Flow Fume Hood (Design Opening)



Frontier Acela_® High Performance Low Flow Fume Hood (Full Close)

Standards Compliance

Chemical Fume Containment

ANSI/ASHRAE 110-1995, USA EN 14175, Europe CAN/CSA-22.2, No.61010-1 EN-61010-1, Europe IEC-61010-1, Worldwide

Electrical Safety







Safety Certified and Tested

All Esco fume hoods are manufactured for the most demanding laboratory applications.

- Independently tested and certified (4ft, 5ft and 6ft models) by Tintschl Engineering AG to the European Standard EN 14175-3.
- Independently tested and certified by Exposure Control Technologies, Inc. to the American Standard ASHRAE 110.
- UL1805 Listed
- Ergonomic design features combine to create a more comfortable work environment, which promotes safety by enabling the user to concentrate on his / her work.
- Hoods are ASHRAE tested on the production line on a sampling basis (industry exclusive).
- Sash is clearly labeled with all operating instructions and illustrations.

8 Superstructure

- · Tri-wall construction for maximum robustness.
- Fabricated of electro-galvanized steel (zinc coated to prevent



rusting even if the powder coat is abraded, and phosphate for better adhesion of the powder coating).

- Epoxy-polyester hybrid Isocide™ antimicrobial powder coating provides long term chemical, abrasion and weathering resistance.
- The pan type design from a single piece of steel (unlike frame welded designs) is easier to assemble on site, and more robust for longterm usage.

- Ergonomic design on the Esco Frontier Acela_® is anthropometrically proven.
- Hood lighting is pre-wired. Lighting is electronically ballasted, energy efficient, instant start. Typical light intensity on work surface is >1076 lux (>100 foot candles) in zero ambient conditions.

Additional Safety Features

- Sash stop limits sash movement beyond 457 mm (18.0"), 'encouraging' user to work at safe positions.
- When sash is raised above 457 mm (18.0"), it will automatically and gently fall back to the safe level (Esco's creep-down mechanism) unless held in place. This enables the sash to be raised temporarily to the full open position for set up of equipment and apparatus in hood, while enforcing




- regular operation of the hood with a lowered sash.
- Sash creep down can only be disabled with a key, which permits
 the laboratory manager or safety officer to restrict operation
 of the hood at sash openings above 457 mm (18.0").
- The sash can also be locked down in the fully closed position with the key, in the event of a fume hood fault, facility shutdown, or, simply to restrict access to equipment set up in the hood.
- Fluorescent light casing is not secured, acting as an explosion relief, allowing upward pressure release during explosions, maximizing user safety.

Aesthetics and Ergonomics Design

- Perfect Pitch Profile™. The gentle pitch of the hood enhances visibility and enables users to work further into the hood without strain.
- A tall sight line creates a bright, airy work space.
- Coordinated color scheme (neutral white with light grey accents) blends the hood with any casework, metal or wood.
- Color temperature is tuned to provide a gentle, bright (but not harsh) and comfortable work environment for the user.

General Specifications, Frontier Acela _® Fume Hood						
Model			EFA-4UDR_W	EFA-5UDR_W	EFA-6UDR_W	EFA-8UDR_W
Nominal Size	е		1.2 meters (4')	1.5 meters (5')	1.8 meters (6')	2.4 meters (8')
External Din (W x D x H)	nensions		1220 x 900 x 1521 mm 48.0" x 35.4"x 59.9"	1525 x 900 x 1521 mm 60.0" x 35.4"x 59.9"	1830 x 900 x 1521 mm 72.0" x 35.4"x 59.9"	2440 x 900 x 1521 mm 96.0" x 35.4"x 59.9"
Internal Dimensions * (W x D x H)						2210 x 675 x 1230 mm 87.0" x 26.6"x 48.4"
Maximum Sash Opening				740 mm	1 (29.1")	
	Face Velocity	Sash Opening				
	0.3 m/s (60 fpm)	457 mm (18.0")	541 cmh @ 14.3 Pa 316 cfm @ 0.06" WG	777 cmh @ 14.6 Pa 457 cfm @ 0.06" WG	872 cmh @ 19.9 Pa 510 cfm @ 0.08" WG	1203 cmh @ 14.2 Pa 708 cfm @ 0.06" WG
Exhaust	0.4 m/s (80 fpm)	457 mm (18.0")	721cmh @ 19.7 Pa 424 cfm @ 0.08" WG	942 cmh @ 23.7 Pa 554 cfm @ 0.10" WG	1163 cmh @ 28.8 Pa 684 cfm @ 0.12 " WG	1604 cmh @ 26.4 Pa 944 cfm @ 0.11 " WG
Volume/ Static	0.5 m/s (100 fpm)	457 mm (18.0")	901 cmh @ 31.8 Pa 530 cfm @ 0.13" WG	1177 cmh @ 34.7 Pa 692.8 cfm @ 0.14" WG	1453 cmh @ 41.8 Pa 855 cfm @ 0.17" WG	2005 cmh @ 32.3 Pa 1180 cfm @ 0.13" WG
Pressure Required	0.3 m/s (60 fpm)	Full	899 cmh @ 22.1 Pa 526 cfm @ 0.09" WG	1175 cmh @ 28.7 Pa 691 cfm @ 0.12" WG	1450 cmh @ 36.1Pa 848 cfm @ 0.15 " WG	1819 cmh @ 27.3 Pa 1070 cfm @ 0.11" WG
	0.4 m/s (80 fpm)	Full	1199 cmh @ 36.7 Pa 701 cfm @ 0.15 " WG	1556 cmh @ 49.3 Pa 922 cfm @ 0.20" WG	1933 cmh @ 61.4 Pa 1138 cfm @ 0.25" WG	2668 cmh @ 48.3 Pa 1570 cfm @ 0.19" WG
	0.5 m/s (100 fpm)	Full	1499 cmh @ 66.6 Pa 877 cfm @ 0.27" WG	1958 cmh @ 76.6 Pa 1152 cfm @ 0.31" WG	2197 cmh @ 94.7 Pa 1285 cfm @ 0.38" WG	3335 cmh @ 74.3 Pa 1962 cfm @ 0.30" WG
Exhaust Out	let Diameter		305 mm (12.0")			
Number of E	xhaust Collars		1	1	1	2

^{*} When sash is in full open position, total height of fume hood is 1644mm/ 64.7". Please ensure sufficient ceiling height clearance is available.



^{*} Combination Sash not available for EFA-8UDR_W-_ model.




Serviceability

- Baffles remove to allow cleaning inside the hood.
- Internal access panels on side walls facilitate access to plumbing connections for service, especially when hoods are installed next to each other.
- Removable front panel facilitates easy access to lighting and other electrical components mounted above the hood work chamber.
- Chain and sprocket sash system requires minimal service.

Warranty

The Frontier $Acela_{\otimes}$ High Performance Fume Hood is warranted for 1 year excluding consumable parts and accessories. Contact your local sales representative for specific warranty details.






maximum opening of the face, but allow


access to the top interior of the hood.

Model EFA, Frontier Acela Fume Hood Technical Specifications

Top View
(for EFA-4UDR_W-_, EFA-5UDR_W-_,
EFA-6UDR_W-_)


Top View
(for EFA-8UDR_W-_)





Front View

10



Side View

(3.6.1.)

900 mm (35.4.")

936 mm (36.9.")

- 1. Fiberglass exhaust collar
- 2. Electrical junction box
- 3. Light housing
- 4. Removable baffle
- 5. Side access panel
- 6. Swan-neck water faucet
- 7. Flash airfoil, hinged for cleaning, spillage containment 0.4 Litres/0.1 gallons
- 8. Removable front panel
- 9. Tempered laminated framed sash glass

- 10. Provision for airflow alarm hinged airfoil
- 11. Light and fan switch
- 12. Distillation grid provision
- 13. Side vane
- 14. Electrical outlets
- 15 Sach ston
- 16. Sash keylock (lock sash at fully open and fully close)
- 17. Provision for additional outlets
- 18. Gas fixture remote
- 19. Cold water fixture remote

Fume Hood Dimensional Data					
Model	A	В	С		
EFA-4UDR_W	1220 mm (48.0")	996 mm (39.0")	610 mm (24.0")		
EFA-5UDR_W	1525 mm (60.0")	1301 mm (51.0")	763 mm (30.0")		
EFA-6UDR_W	1830 mm (72.0")	1606 mm (63.0")	915 mm (36.0")		
EFA-8UDR_W	2440 mm (96.0")	2210 mm (87.0")	664 mm (26.0")		

Fume Hood Installation Requirements

- Proper location (refer to Esco recommendations)
- Exhaust system capable of delivering airflow and pressure drop requirements
 - Exhaust ductwork and connection accesories
 - Exhaust blower (or building exhaust system)
 - Contactor (if applicable)
 - Means of adjusting exhaust flow (damper or frequency inverter)
- Sufficient laboratory supply air
- Base cabinet or other means of support
- One set of filler panels for each continuous row of hoods
- Plumbing and electrical hook-up
- Unless otherwise specified, items and services specified above are not provided by Esco


Frontier, **acela**™

Options and Accessories

Esco offers a variety of options and accessories to meet application requirements.


All Frontier Acela® Fume Hoods include 1 water and 1 gas Enhanz™ remote-controlled fixture plumbed to the top of hood and and 4 single outlets standard. (Specify outlet code when ordering)





Enhanz™ **Service Fixtures**

Remote control service fixtures for fume hoods, preplumbed. Up to 8 service fixtures per fume hood.



A: Argon **HW**: Hot Water CO2 : Carbon Dioxide DI: Deionised Water

Sentinel•XL[™] **Airflow Alarm Kit**

Real-time digital display of face velocity. Audible and visual alarms.

SXL – EFA



Distillation Grids

Scaffoldings to support clamps for distillation apparatus.

DG - 4U S

External Width

Material 4U: 1220mm (48") **S**: Stainless steel

5U: 1525mm (60") **6U**: 1830mm (72") **8U**: 2440mm (96")

Drip-Cups


SPP - DC

Material

SPP: Polypropylene SIC : Ceramic,

Only For Ceramic Work Tops





Туре	SS304	SS316	Trespa TopLab ^{Plus} Phenolic Resin	Ероху	Ceramic
Chemical Resistance	Good	Better	Better	Better	Best
Temperature Resistance	Better (300°C)	Better (300°C)	Good (110°C)	Good (165°C)	Best (1200°C)
Cost	Low Price	Mid Price	Mid Price	Mid Price	Premium

Fume Hood Work Surface Selection Esco offers 5 types of fume hood work surfaces for different applications



Frontier Acela_® Base Cabinetry (EBA)



Safety Certified and Tested

 Built and tested according to SEFA-8 recommended practices.

Superstructure

 Fabricated of electro-galvanized steel (zinc coated to prevent rusting even if the powder coat is abraded, and phosphated for better adhesion of the powder coating).

- Epoxy-polyester hybrid Isocide™ antimicrobial powder coating provides long term chemical, abrasion and weathering resistance.
- Ships unassembled, assembles rapidly on site.

Aesthetics and Ergonomics Design

- Concealed door hinges.
- Soft-close insulated doors.
- Adjustable shelf allows user to select optimum heights for upper and lower storage compartments.

Serviceability

- Convenient access to levelers from inside base cabinet.
- Removable rear access panel facilitates installation and maintenance of plumbing and drainage systems.

Options and Accesories

MCB/ELCB

- Protects laboratory equipment during sudden fluctuation of current.
- Fume hood circuit protection.
- Only applicable to 230V, AC, 50/60Hz hoods.
- Factory-installed; specify when ordering.

Ventilation Kit (VK-EBA)

- Ventilates base cabinet utilizing the hood exhaust system.
- Field-installed.

• Filler Panel (FP-EBAD)

- One set of filler panels required per continous row of hoods.
- Field-installed






No.	Type of Test	Test Result
1	Cabinet load test	PASS
2	Cabinet concentrated load test	PASS
3	Cabinet torsion	PASS
4	Cabinet submersion test	PASS
5	Door hinge test	PASS
6	Door impact test	PASS
7	Door cycle test	PASS
8	Chemical spot test	PASS
9	Hot water test	PASS
10	Impact test	PASS
11	Paint adhesion on steel	PASS
12	Paint hardness on steel	PASS



Door Hinge Test


Cabinet Load Test

Model EBA, Frontier Acela_® Base Cabinet Technical Specification




- 1. Rear access panel
- 2. Soft close door
- 3. Adjustable shelf
- 4. Adjustable shelf slot
- 5. Base cabinet height leveler
- 6. Base cabinet filler panel (optional)
- 7. MCB/ELCB set (optional)
- 8. Ventilation provision

Front View EBA-4UDG-_ (for EFA-4UDR_W-_)




Front View


EBA-5UDG-_/EBA-6UDG-_/EBA-8UDG-_
(for EFA-5UDR_W-_, EFA-6UDR_W-_, EFA-8UDR_W-_)



Top View (for All Sizes)



Side View (for All Sizes)



Model*	A	В
EBA-4UDG-0 /	1220 mm	610 mm
EBA-4UDG-8	(48.0")	(24.0")
EBA-5UDG-0 /	1525 mm	381 mm
EBA-5UDG-8	(60.0")	(15.0")
EBA-6UDG-0 /	1830 mm	457.5 mm
EBA-6UDG-8	(72.0")	(18.0")
EBA-8UDG-0 /	2440 mm	610 mm
EBA-8UDG-8	(96.0")	(24.0")

* Model EBA-____-0 is standard base cabinet without MCB/ELCB.

Model EBA-____-8 is with MCB/ELCB protection and only applicable for 230V,AC,50/60Hz units.

MCB/ELCB panel is always mounted on the top right hand corner of base cabinet set.



Esco Resinate™

Esco Resinate™ is a proprietary composite material specifically designed for use as internal liner in laboratory fume hoods.

- Excellent chemical resistance (refer to table below)
- Excellent physical properties provide structural reinforcement for the hood
- Smooth, attractive, easy-to-clean finish


	Chemic	cal Resistance* of Es	sco Resinate™	Internal Liner		
	Chemicals	Result	_	Chemicals	Result	
	85% Sulfuric Acid	No Effect		40% Sodium Hydroxide	No Effec	
	98% Sulfuric Acid	1st Grade	Alkalis	65% KOH	No Effec	
	50% Nitric Acid	1st Grade		10% Iron Chloride	No Effec	
	65% Nitric Acid	2nd Grade		10% Copper Sulfate	No Effec	
Acids	36% Hydrochloric Acid	No Effect		15% Sodium Sulfide	No Effec	
	85% Phosphoric Acid	No Effect		Ammonium Hydroxide	No Effec	
	40% Hydrofluoric Acid	No Effect	_			
	60% Chromic Trioxide	No Effect		Chemicals	Result	
	99% Glacial Acetic Acid	No Effect	General Reagents	50% Magnesium Sulfate	No Effec	
	Aqua Regia	No Effect		34% Hydrogen Peroxide	No Effec	
_				Urea	No Effec	
	Chemicals	Result		Copper Sulfate	No Effec	
	37% Fomaldehyde	No Effect		Karl Fisher Reagent	No Effec	
	N-Hexane	No Effect		Iodine	No Effec	
	Ethyl Acetate	No Effect	_			
	Ethyl Ether	No Effect		Chemicals	Result	
	Ethyl Alcohol	No Effect		1%Gentian Violet	No Effec	
	Isopropyl Alcohol	No Effect	Stains and	Methylene Blue	No Effec	
olvents	Carbon Tetrachloride	No Effect	Indicators	Crystal Violet	No Effec	
orvents	Naphthalene	No Effect		Methyl Red	No Effec	
	Chloroform	No Effect		Methyl Orange	No Effec	
	Methanol	No Effect	_			
	Toluene	No Effect	Note: Esco Resinate™ may not be suitable for fume hoods for increased acidic and heat loads. Contact your local sales representative for details. * Test Method: One drop of test chemical placed on material surface and covered			
	Xylene	No Effect				
	Acetone	No Effect				
	Styrene	No Effect	with watch gl	ass for 16 hours before result is observed.		
	Phenol	No Effect	 1st Grade: Slight effect on color and gloss. No change to physical properties. 2nd Grade: Clear effect on color and gloss. No change to physical properties. 			

Esco Resinate PlusTM

Esco Resinate Plus liner is offer excellent chemical and physical resistance against harsh environments particularly against highly corrosive acids.

- Fiberglass Reinforced Plastic
- UL1805 Compliant
- Smooth, attractive, easy-to-clean finish

Frontier, ACELA"







Since 1978, Esco has emerged as a leader in the development of controlled environment, laboratory and pharmaceutical equipment solutions. Products sold in more than 100 countries include biological safety cabinets, compounding pharmacy equipment, containment I pharma products, ductless fume hoods, in vitro fertilization workstations, lab animal research products, laboratory fume hoods, laboratory ovens and incubators, laminar flow clean benches and PCR products and instrumentation. With the most extensive product line in the industry, Esco has passed more tests, in more languages, for more certifications, throughout more countries than any biosafety cabinet manufacturer in the world. Esco remains dedicated to delivering innovative solutions for the clinical, life science, research and industrial laboratory community. www.escoglobal.com.

Life Science • Chemical Research • Assisted Reproductive Technology (ART) • Pharmaceutical Equipment • General Equipment



WORLD CLASS. WORLDWIDE.

Esco Technologies, Inc. • 2940 Turnpike Drive, Units 15-16 • Hatboro, PA 19040, USA Toll-Free USA and Canada 877-479-3726 • Tel 215-441-9661 • Fax 215-441-9660 us.escoglobal.com • usa@escoglobal.com

Esco Micro Pte. Ltd. • 21 Changi South Street 1 • Singapore 486 777 Tel +65 6542 0833 • Fax +65 6542 6920 • mail@escoglobal.com www.escoglobal.com

Esco Global Offices | Beijing, China | Chengdu, China | Kuala Lumpur, Malaysia | Manama, Bahrain | Guangzhou, China | Hanoi, Vietnam | Melaka, Malaysia | Mumbai, India | Philadelphia, PA, USA | Salisbury, UK | Shanghai, China | Seoul, Korea | Delhi, India | Osaka, Japan | Manila, Philippines | Midrand, South Africa | Jakarta, Indonesia | Singapore





