

# Releasat®

## BIOLOGICAL INDICATOR CULTURING SET

*For Industrial Use Only*

### CERTIFICATE OF ANALYSIS

Reorder No: RCD/50

*Bacillus atrophaeus* 9372<sup>(1)</sup>

Biological Indicator for: Chlorine Dioxide sterilization

Culture: 36 - 38°C. The supplied bacteriological medium will meet requirements for growth promoting ability.

Purity: No evidence of contaminants using standard plate count techniques.

Releasat Gas Lot No.: **RCD-**

Spore Strip Lot No: ACDS-

Media Lot No: CM-

Manufacture Date:

Expiration Date: 12 months from Manufacture Date

Heat Shocked Population: x 10<sup>6</sup> Spores / Unit

Carrier Size: ¼" x ¾" (6 mm x 19 mm)

Units are manufactured in compliance with Mesa Laboratories' quality standards, USP, and ISO 11138 guidelines and all appropriate subsections.

<sup>(1)</sup>Culture is traceable to a recognized culture collection identified in USP and ISO 11138.

Certified By: \_\_\_\_\_  
Quality Representative



**Bozeman Manufacturing Facility**

10 Evergreen Drive

Bozeman, MT 59715

T: 406/585-9535 ♦ F: 406/585-9219

www.sgmbiotech.com ♦ www.mesalabs.com